



MSF Newborn Warming Project Market Survey: Warming Devices

Saige Sunier, Timothy Prestero - Design that Matters

Overall Summary

[illegible]

Product Requirements: Portability and CE Certification

Conductive Warming Mattresses

Ardo Amecosy



Mediprema Calibed



Medwarm IM-55M



AstroPad Warming Blanket



Kanmed BW3-001



Inspiration CosyTherm²



Thermocare Convenience



Ginevri Aquatherm



Ramonak-01



Chemical-Thermal Batteries

Warmilu IncuBlanket



Warmilu Wipeable IncuBlanket



Embrace Nest Warmer



DreamWarmer



TransWarmer



BabyPod II
(with TransWarmer)



Other

Ceratherm Radiant Warmer




Transport Incubator*



*Portability requirement implies hand-carry rather than on wheels, given the need to transport patients up stairs, through narrow doorways and crowded wards, and outside over rough ground

Product Requirement: Context Appropriate (Durable and Rugged, Limited Maintenance)

Failure Mode (What/When)



Fragile cable between heating source and control unit are prone to damage.

Evidence (Why/Who)

Goyalmara Field Staff: The cable between the Amecosy mattress and the control unit breaks easily and cannot be repaired or replaced.

Yemen Pediatric Staff: It is very difficult to get Biomed Staff into Yemen, it is important that the device requires minimum and easy maintenance.

Stakeholder Survey: 75% rank as a "Must Have"; 25% rank as a "Nice to Have"

Solutions (How)

1. Control unit and mattress are one unit with no connecting cords.
2. Connector cable is easily replaced by locally available consumables.
3. Product is installed in a bassinet or in one area and rarely moved.

Devices Affected

Ardo Amecosy



Mediprema Calibed



Medwarm IM-55M



AstroPad Warming Blanket



Warmilu IncuBlanket



Warmilu Wipeable IncuBlanket



Embrace Nest Warmer



Ceratherm Radiant Warmer



Ramonak-01



CosyTherm²



Kanmed BW3-001



Thermocare Convenience



Ginevri Aquatherm



DreamWarmer



TransWarmer



BabyPod II (with TransWarmer)



Transport Incubator



Product Requirement: Context Appropriate (IPC: Easy to clean between patients, with common cleaning products)

Failure Mode (What/When)



Difficult and/or extensive process to sterilize and clean product between patients. Is not compatible with common cleaning products.

Evidence (Why/Who)

Workshop Feedback: Some sites have a very high patient to staff ratio and the patient turnover is very rapid, making a quick easy IPC process important.

Workshop Feedback: Sites will use bleach and chlorine products to sterilize equipment.

Stakeholder Survey: 100% rank as "Must Have"

Solutions (How)

1. Choose a product that can be wiped down for sterilization between patients.
2. Product must have minimal creases and seams to reduce the risk of collecting dirt and infectious germs.

Devices Affected

Ardo Amecosy



Mediprema Calibed



Medwarm IM-55M



AstroPad Warming Blanket



Warmilu IncuBlanket



Warmilu Wipeable IncuBlanket



Embrace Nest Warmer



Ceratherm Radiant Warmer



Ramonak-01



CosyTherm²



Kanmed BW3-001



Thermocare Convenience



Ginevri Aquatherm



DreamWarmer



TransWarmer




BabyPod II (with TransWarmer)



Transport Incubator




Product Requirements: (1) Context Appropriate (Limited infrastructure reqs); (2) User Friendly (no risk of scalding or hyperthermia)

Failure Mode (What/When)	Evidence (Why/Who)	Solutions (How)
<div><p>Requires room temperature control (AC) for accurate and safe performance.</p></div>	<p>Yemen Pediatric Staff: The project site can exceed 40C.</p> <p>Workshop Feedback: Devices in the catalog will be used in all MSF sites and should be usable in all climates.</p> <p>Stakeholder Survey (1): 75% rank as a "Must Have"; 25% rank as a "Nice to Have"</p> <p>Stakeholder Survey (2): 100% rank as "Must Have"</p>	<ol style="list-style-type: none">1. If the product is a warming mattress the temperature probe must read from the surface of the mattress and record real time temperature.2. Operating the warmer (whether chemical-thermal battery or electrically-powered device) at an ambient temp >28C does not introduce a risk of scalding.

Devices Affected									
									
									

Product Requirements: Context Appropriate (Low Cost of Ownership, Limited Consumables)

Failure Mode (What/When)



High consumable requirement cannot be maintained in all MSF project sites.

Evidence (Why/Who)

Yemen Pediatric Staff: It is very difficult to get Biomed staff, medical materials, and medical products into Yemen because of the political situation.

Kusisa, DRC Staff: The project site is very difficult to reach and any consumables will need to be hand carried or carried on a motorcycle. The political situation makes transport risky.

Working Group: While it is not MSF's priority, any device that relies on consumables is not sustainable once MSF leaves.

Stakeholder Survey: 25% rank as "Must Have"; 63% rank as "Nice to Have"; rank as 13% "Don't Know"


Solutions (How)

1. Choose a product that requires very little, if any, consumables.
2. If product requires consumables, the parts can be purchased locally.
3. Provide additional method of warming if choosing a consumable warming device.

Devices Affected



Product Requirement: (1) User Friendly (Device state is obvious); (2) User Friendly (limited risk of scalding or hyperthermia)

Failure Mode (What/When)	Evidence (Why/Who)	Solutions (How)
<div><p>Device state (on, off, broken) is not obvious. Device can break or temperature become too high/low without user knowing.</p></div>	<p>Working Group: Some MSF project sites have a high patient: staff ratio which makes it difficult to check device status when focusing on patient care.</p> <p>Stakeholder Survey (1): 75% rank as “Must Have”; 25% rank as “Nice to Have”</p> <p>Stakeholder Survey (2): 100% rank as “Must Have”</p>	<ol style="list-style-type: none">1. Device has auditory and visual alarms that notify user if device is not functioning as expected.2. Temperature feedback control: device automatically monitors and adjusts temperature to maintain the appropriate range, or shuts off if ambient temperature is too high.

Devices Affected

